

BTEX

Primary reference materials containing **BTEX in nitrogen** in specific concentrations. Within the listed ranges any composition can be supplied.

PRM code	Certified component	Nominal fraction range (mol/mol)	Stability period
BTEX	Benzene, Toluene, Ethylbenzene, o-Xylene, m-Xylene, p-Xylene	10 - 100 · 10 ⁻⁹	24 months

Cylinder pressure 11,5 MPa, Cylinder valve outlet DIN 477 no. 1, Approximate deliverable volume 0,6 m³.