

# Carbonyl sulphide in methane

Primary reference materials containing **carbonyl sulfide in methane** are available in a broad range of composition. Within the listed ranges any composition can be supplied.

PRM code	Certified component	Nominal fraction range (mol/mol)	Stability period
JC21	COS	25 - 250 · 10 <sup>-6</sup>	24 months
JC22	COS	250 - 2500 · 10 <sup>-6</sup>	24 months

Cylinder pressure 11,5 MPa, Cylinder valve outlet DIN 477 no. 1, Approximate deliverable volume 0,6 m<sup>3</sup>.