

Nitric oxide

Primary reference materials containing **nitric oxide in nitrogen** are available in a broad range of composition. Within the listed ranges any composition can be supplied. Other concentrations and other matrices are available on request.

PRM code	Certified component	Nominal fraction range (mol/mol)	Stability period
BD09	NO	0,1 - 1 · 10 ⁻⁶	24 months
BD10	NO	1 - 10 · 10 ⁻⁶	24 months
BD11	NO	10 - 100 · 10 ⁻⁶	36 months
BD12	NO	100 - 1000 · 10 ⁻⁶	36 months
BD13	NO	0,1 - 0,8 · 10 ⁻²	36 months

Cylinder pressure 11,5 MPa, Cylinder valve outlet DIN 477 no. 1, Approximate deliverable volume 0,6 m³.