

# Nitrous oxide

Primary reference materials containing **nitrous oxide in synthetic air** in specific concentrations.

Other concentrations and other matrices are available on request.

PRM code	Certified component	Nominal fraction range (mol/mol)	Stability period
BI21	N <sub>2</sub> O	0,3 - 3,0 · 10 <sup>-6</sup>	24 months

Cylinder pressure 11,5 MPa, Cylinder valve outlet DIN 477 no. 1, Approximate deliverable volume 0,6 m<sup>3</sup>.