

Water

Primary reference materials containing **moisture in nitrogen** are available in the listed range.

Within the listed range any composition can be supplied.

| PRM code | Certified component | Nominal fraction range (mol/mol) | Stability period |
|----------|---------------------|----------------------------------|------------------|
| PRM code | Certified component | Nominal fraction range (mol/mol) | Stability period |
| BC11 | H ₂ O | 25 - 200 · 10 ⁻⁶ | 24 months |

Cylinder pressure 7,0 MPa, Cylinder valve outlet DIN 477 no. 1, Approximate deliverable volume 0,4 m³.